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For research or manufacturing use only.

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5PRP Kit V1 kit is used to isolate and activate platelet-rich plasma from 12ml of peripheral blood, umbilical cord blood, etc.

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5PRP Kit V1 is a kit for blood collection and manufacturing of activated platelet-rich plasma. This kit draws blood using vacuum tubes and separates platelets by isopycnic centrifugation. PRP is activated through a filtration membrane, resulting in a final activated PRP product that does not contain any blood cells and has a high purity level. The product manufacturing process is done quickly with 1 single centrifugal step. 5PRP Kit V1 was researched and developed by Stem Cell Institute,

HCMUS, VNUHCM, Viet Nam. 5PRP Kit V1 guarantees to be sterile (0 CFU) with endotoxin level ≤ 1 EU/mL..

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5PRP Kit V1 provides a fast and convenient method to collect blood and manufacture platelet-rich plasma from peripheral blood or umbilical cord blood, etc.

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Storage and transportation temperature: Room temperature.

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Expiry date: following the nearest expiry date of the components.

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Attach the Vacuum blood collection needle (butterfly needle) to the Needle holder.

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Clean the puncture site with an alcohol pad, then perform venipuncture.

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Advance the PRP tube into the Needle holder; at this point, blood will flow automatically into the tube.

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Collect blood until the tube is full (tube contains 10-12 mL of blood). Gently invert tubes to mix blood with an additive.

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Place the PRP tubes containing blood into the centrifuge (place 1 more counterweight tube in a symmetrical position).

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Centrifuge the tubes at 3500 RPM for 10 minutes.

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Gently take out the PRP tubes from the centrifuge.

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Attach the 10 mL Luer-lock Syringe to a Spinal anesthesia needle. Remove supernatant on the top and leave 3-5 mL of plasma at the bottom (containing platelet pellet).

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Suspend the platelet pellet in 3-5 mL of plasma, and aspirate 3-5 times. This mixture is inactivated platelet-rich plasma (PRP).

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Collect inactivated PRP into the 10 mL Luer-lock Syringe.

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Attach the 10 mL Luer-lock Syringe to the Filter.

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Attach the 5 mL Luer-lock Syringe to the other end of the Connector.

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Strongly push the piston of the 10 mL Luer-lock Syringe so that inactivated PRP goes through the Filter and into the 5 mL Luer-lock Syringe (The product in 5 mL Luer-lock Syringe is activated PRP).

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Do not use the kit if the packaging is damaged or cracked, or if one of the kit components is abnormal.

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The symbols present on the kit are explained below

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